

# Command:SPIN (@)

MS2000 or RM2000 syntax

<b>Shortcut</b>	@
<b>Format</b>	SPIN [axis]=[-128 to 128]...
<b>Units</b>	Integer

Tiger syntax

<b>Shortcut</b>	@
<b>Format</b>	SPIN [axis]=[-128 to 128]...
<b>Units</b>	Integer
<b>Type</b>	Axis-Specific

Tells controller to 'spin' the motor of specified axis at a rate expressed as its DAC value, a bit value from -128 to 128. This causes the motor to run in open-loop mode with no position or speed feedback.

## Reply

If there are no errors, a positive reply of :A is sent back.

## Example

```
@ X=100 Y=-100 Z  
:A
```

This example shows a command that will instruct the X-axis turn at a motor rate of 100 DAC bits in one direction, the Y-axis at the same rate but in the other direction, and stop any rotation or motion of the Z-axis.

**Note 1:** To stop rotation, give a value of zero, or just the type the axis letter without an assignment as shown in the example above, or use the [HALT command](#).

**Note 2:** The HALT command will not return an :N-21 when stopping a SPIN command.

[commands](#), [tiger](#), [ms2000](#), [SPIN](#), [VECTOR](#)

From:

<http://asiimaging.com/docs/> - **Applied Scientific Instrumentation**

Permanent link:

<http://asiimaging.com/docs/commands/spin>

Last update: **2025/04/08 17:59**

