

# Atmospheric Specifications for ASI Tiger Controllers

<b>Operating</b>	
<b>Temperature</b>	10° to 35° celsius
<b>Humidity</b>	0-75% RH non condensing
<b>Atmospheric Pressure</b>	700 HPa to 1060 HPa
<b>Storage</b>	
<b>Temperature</b>	-20° to 45° celsius †
<b>Humidity</b>	10-90% RH non condensing
<b>Atmospheric Pressure</b>	500 HPa to 1060 HPa

(†) must let product return to room temperature over 24 hour period before operating.

**Note 1:** The above specs are for the ASI Tiger Controllers. Other ASI parts may have different Specifications - please inquire for individual part specs.

**Note 2:** These are calculated specs going on ratings from individual components utilized. ASI does not test extreme conditions.

[safety](#)

From:  
<http://asiimaging.com/docs/> - **Applied Scientific Instrumentation**

Permanent link:  
[http://asiimaging.com/docs/atmospheric\\_spec](http://asiimaging.com/docs/atmospheric_spec)

Last update: **2026/06/05 03:59**

