

FTP Controller



The toggle switch near the joystick provides for a rapid means of moving the stage up & down. Fine control of z position is accomplished with the Z focus knob, or by computer control.



There is an option to control the Z focus with foot switches as well. If this option is requested the 9 pin serial out connector is used to connect the foot switches, and is labeled as such.

[FTP Setup Manual](#)

[FTP Firmware Manual](#)

[manual, ms2000, ftp](#)

From: <https://asiimaging.com/docs/> - **Applied Scientific Instrumentation**

Permanent link: https://asiimaging.com/docs/ftp_controller

Last update: **2026/06/19 17:21**

