


Olympus IX71 Linear Encoder Installation



Please excuse the editing. This document is still being drafted.



Linear encoder holder mounts via single 5×30 mm screw.



Space encoder to be almost fully depressed when at focal plane. This will allow maximum travel in the reverse direction.

[olympus](#), [ix71](#), [linear encoder](#), [zdrive](#)

From:

<https://asiimaging.com/docs/> - **Applied Scientific Instrumentation**

Permanent link:

https://asiimaging.com/docs/olympus_ix71_linear_encoder_install

Last update: **2021/09/23 17:15**

