

# Solenoid Lock Insert (I-3020-MLLAC)

I-3020-MLLAC is a solenoid powered insert that is able to secure a Multi-well plate in place.

The insert is driven with ASI MS2000 or TG-1000 controller.



Solenoid Lock Insert (I-3020-MLLAC)



Insert's connector a 4pin HiRose connector



Corresponding connector on a MS2000 controller

The insert is operated with the [SECURE command](#)  
[ms2000](#), [tiger](#), [manual](#), [solenoid lock](#)

From:

<http://asiimaging.com/docs/> - **Applied Scientific Instrumentation**

Permanent link:

[http://asiimaging.com/docs/solenoid\\_lock\\_insert](http://asiimaging.com/docs/solenoid_lock_insert)

Last update: **2021/09/23 17:15**

